

Visagraph Installation Manual


Visagraph Installation Apple (Read Only)

IMPORTANT: The Apple version of Visagraph is read only. At this time it cannot be used to take recordings.


Go to <https://visagraph.readingplus.com/> to download the Visagraph Application.

If you have any questions regarding the installation of Visagraph, contact the Reading Plus Customer Service department. Our contact information can be found [here](#).


Click “Finder”




Click on “Downloads”




Double Click “Visagraph_Installer_Mac_Read-Only_V1.0.dmg”




Drag Visagraph (Read Only) on top of the Applications Icon.



Drag Visagraph (Read Only) to your desktop.




Double Click the Visagraph Icon



Enter your username and password.

(Optional) – Check the box “Remember username and password”

Click “Login”

A screenshot of a web browser window titled "Visagraph (Read Only)". The browser's address bar shows the Visagraph logo and "Taylor Associates". The main content area has a light blue background. In the center, the text "Visagraph Login" is displayed in a large, dark blue font. Below this text are two horizontal yellow input fields for username and password. Under the password field is a checkbox labeled "Remember username and password". At the bottom center is a white button with the text "Login".